

# VIRAL PARTICLE PURIFICATION

## Reagents

- Light CsCl (1.2 gm/mL in 10 mM tris pH 8.1)
- Heavy CsCl (1.45 gm/mL in 10 mM Tris pH 8.1)
- 10 mM Tris pH 8.1
- 3% Sucrose/PBS (3 L for Dialysis per virus)
- 50% Glycerol/BSA

## Materials

- Slide-A-Lyzer Cassettes from Pierce – Cat #:66425
- Beckman Centrifuge Tubes, 1 x 3 inches – Cat #:344058
- Beckman Centrifuge Tubes,  $\frac{9}{16}$  x 3 inches – Cat #:344059

## Procedure

- Freeze/thaw the lysate 3 times, using a dry ice EtOH bath and 37°C waterbath.
- Spin the lysate at 3500 rpm ( $g = 1177.8$ ) for about 10 min.
- First Cesium Chloride gradient.
  - a. 20 plate lysate -
    - ◆ Sw 28 rotor and buckets for Beckman Optima LE-80K ultracentrifuge.
    - ◆ Add 12 ml of light CsCl to Beckman centrifuge tube.
    - ◆ Add 12 ml of heavy CsCl very slowly under the light CsCl. Take 14 ml of heavy CsCl in the pipette so no air bubbles disrupt the gradient.
    - ◆ Add the clarified lysate very gently against the side of the tube on top of the CsCl gradient. Supplement with 10 mM Tris to about  $\frac{1}{4}$ " below the top of the tube. Make sure that the volume is balanced between the 2<sup>nd</sup> sample and/or the blank.
  - b. 10 plate and 5 plate lysate -
    - ◆ For Beckman Ultras, use the Sw 41 Ti set.
    - ◆ Add 3 ml of light CsCl to the centrifuge tube.
    - ◆ Add 3 ml of heavy CsCl very slowly under the light CsCl. Always take 4-6 ml of heavy CsCl in the pipette so no air bubbles disrupt the gradient.

- ◆ Add the clarified supernatant very gently against the side of the tube on top of the CsCl gradient. Supplement with 10 mM Tris to about  $\frac{1}{4}$ " below the top of the tube. Make sure the tubes are balanced.
- Spin in the ultra centrifuge at 20K rpm ( $g = 33,735$ ) at  $4^{\circ}\text{C}$  for a minimum of 3 h.
- Pull the viral band. The viral band is the bottom white colored band that is present where the heavy and light CsCl meet:
  - a. Remove the cell debris at the top of the gradient, to about half an inch above the viral bands with a glass Pasteur pipette.
  - b. Place solution in waste bottle.
  - c. Use a 5 cc syringe with an 18 g 1.0" needle. Insert the needle through the side of the tube about four needle widths below the virus band.
  - d. Extract the virus in as small a volume as possible, leaving behind incomplete bands higher up the tube.
  - e. Dilute the band by raising the volume to 3 ml of 10 mM Tris in a snap cap tube. Place the syringe and needle inside the snap cap tube and pull up the 10 mM Tris. The band must be diluted by 50% with 10 mM Tris.
  - f. If you started with the 5 or 10 plate lysates, dilute your pulled band to only 2 mL.
- Second gradient -
  - a. 20 plate lysate -
    - ◆ Use the Beckman Sw 41 Ti rotor and buckets.
    - ◆ Add 3 ml of light CsCl to centrifuge tube.
    - ◆ Add 3 ml of heavy CsCl very slowly under the light CsCl.
    - ◆ Add the virus/10 mM Tris very gently against the side of the tube on top of the CsCl gradient. Supplement with 10 mM Tris to about  $\frac{1}{4}$ " below the top of the tube.
  - b. 10 or 5 plate lysate -
    - ◆ Use the Beckman Sw 60 or equivalent rotor, buckets and tubes.
    - ◆ For 1 mL gradients, following the above instructions to form the gradient and add the virus.
- Spin in the ultra centrifuge at 20K rpm ( $g = 30,195$ ) at  $4^{\circ}\text{C}$  overnight.
- Pull the viral band as above using a 5 cc syringe. If not dialyzing, put the virus in a small, labeled snap cap tube, or eppendorf microcentrifuge tube.
- Store at  $-80^{\circ}\text{C}$
- Dialysis -
  - a. In the syringe, add 1 ml of 3% Sucrose/PBS, about half the volume of the band. If the virus is too concentrated, it will precipitate out.
  - b. Inject the contents of the syringe into a "Slide-A-Lyzer" dialysis cassette from Pierce. Using the 5 cc syringe pull out the remaining air in the slide.

- c. Place the dialyzer slide with a buoy in a bucket of 3% Sucrose/PBS. Use 1L Sucrose/PBS per slide.
- d. Change the Sucrose/PBS once an hour, 2 times. You will need a total of 3 liters of 3% Sucrose/PBS for each virus prep.
- e. To remove the virus from the dialyzer slide, use a 5 cc syringe and add 1 - 2 mL of air to the slide. Use a different port from step "b".
- f. With syringe still in the Slide-A-Lyzer extract the virus.
- g. Note total volume of virus in the syringe. Subtract 0.4 ml from the total volume.
- h. Place virus in a snap cap tube.
- i. For glycerol stock -
  - ◆ Transfer 250  $\mu$ l of virus into a 2 ml screw cap tube.
  - ◆ Label the tube with virus name.
  - ◆ Add approximately 300  $\mu$ l of 50% Glycerol/BSA to the virus. (For very thick bands, add 500 to 1000  $\mu$ L.)
- j. Turn the spectrophotometer on and allow it to warm up.
- k. Prepare a blank of 13  $\mu$ l of Sucrose/PBS in 987  $\mu$ l of 10 mM Tris.
- l. Dilute the virus to 1:75 dilution (13  $\mu$ l of virus in 987  $\mu$ l of 10 mM Tris).
- m. Measure the absorbance at 260 nm wavelength.
  - ◆  $1.0 \text{ O.D.}_{260} = 1.1 \times 10^{12} \text{ pt/ml}$
  - ◆ Calculate pt/ml concentration taking into consideration the 1:75 dilution factor.
- n. Dilute in 3% Sucrose/PBS or Glycerol/PBS to the virus to  $1.0 \times 10^{12}$  particles/ml.
- o. Target values -
  - ◆ 3% Sucrose/PBS preps:  $1 \times 10^{12}$  pt/ml for the working stock (between 1 - 1.5 is ok).
  - ◆ Glycerol/BSA stocks:  $3 \times 10^{12}$  pt/ml (between 3 - 4 is ok).
- p. Aliquot and freeze at  $-80^{\circ}\text{C}$  (sucrose/BSA) or  $-20^{\circ}\text{C}$  (glycerol/BSA)
  - ◆ For 3% Sucrose/PBS preps:
    - 1 mL, 500  $\mu$ L, 200  $\mu$ l aliquots

- - - o 0 o - - -